

## Table of Contents

1. Overview
2. Workflow
3. Glossary
4. Status Table

## 1. iMITS Overview

iMITS coordinates the production of mutant mice for high-throughput production pipelines to minimize overlap and maximize efficiencies. iMITS captures the intention of a particular consortium to produce mutant mice on a particular gene, the progress of mouse production and the progress of modified allele creation on the original mouse (leading up to the capture of phenotype data on the mouse).

iMITS serves mouse production facilities from around the world as part of its funded mission to support the KOMP2 project- a NIH Common fund project to generate and phenotype 2500 knockout mouse strains in the next five years (1U54HG006370). Mouse production and phenotype centers in the KOMP2 project are required to track mouse production through iMITS to minimize overlap and maximize efficiencies. KOMP2 production is coordinated with other mouse production centers through participation in the International Mouse Phenotyping Consortium (IMPC) whose goal is to generate a knockout mouse strain for every protein-coding gene in the mouse genome in the next 10 years. To ensure a common genetic background, IMPC partners use the targeted ES cells generated by the International Knockout Mouse consortium (IKMC). The majority of IMPC production centers use iMITS to track production to help prevent duplicate strain production and identify how their progress compares to other centers.

The information in iMITS is provided by each individual production center/consortia through the iMITS graphical interface. This is a trust system- iMITS does not have access to a center's laboratory management system and thus cannot double check entries. Therefore it's crucial that personnel entering in the information understand how the interfaces work. In this document we will provide a step-by-step guide. A glossary explains unique terms and a Status Table explains what is meant when iMITS assigns a status to a particular field.

## 2. WORKFLOW

## Gene Selection

## Plans


### Figure 1.

- Workflow for iMits users from registering an interest to a gene to the entering of a mouse strain into a phenotyping pipeline.
- Dashed boxes indicate what iMits module(s) are used to enter in

## Create MI Attempt

## Mouse Production

## Phenotyping



## 2A. Gene Selection: Selecting a gene for Mouse Production.

Users acting on behalf of consortia select genes for mouse production. Users will be alerted if mouse strains with targeted mutations of this gene are in production or have been produced/phenotyped by other IMPC partners. The process works in the following manner:

Users create a plan for the gene(s) of interest. When a plan is created, iMITS will check whether any other consortia have plans (or mice) for the gene. If not, the plan will be immediately put in status “Assigned”. If there are other existing plans or mice, iMITS will put the plan in status “Inspect – Conflict” or “Inspect – MI Attempt” or “Inspect – GLT Mice” depending on the known products.

Conflict statuses do not preclude the centre from actually acting (starting ES-QC or starting mouse production – see below). If the centre does start activity, the plan will be changed by iMITS to have status “Assigned – ES QC in progress” etc.

**NOTE: iMits’ role is not to stop production centres from creating duplicate mouse strains but rather to inform production centres of similar efforts so they may make informed decisions. Duplicate strain production is to be minimized under the KOMP2 funding mechanism.**

1. Select Gene(s) of interest. Here I have selected 1110002L01Rik and filled in the Consortium and Production Centre and Priority drop-downs.

**Register Interest in Micro-Injecting Genes**

Filter Genes

Consortium: BaSH Production Centre: WTSI Sub Project: Priority: High Phenotype Only: Mutagenesis Via: Crispr/Cas9

Register Interest


View In IMPC	Production History	Tree	# KMC Projects	# Clones	Non-Assigned Plans	Assigned Plans	Aborted MIs	MIs In Progress	Genotype Confirmed N	Phenotype Attempts
<a href="#">110001G20Rik</a>	Production...		2	14 Knocko... 6 Targeted... 0 Deletion		<a href="#">[UCD-KOMP.U...</a>	<a href="#">[DTCC:UCD:1]</a>			
<a href="#">1110002J07Rik</a>	Production...		1	4 Knockout... 0 Targeted... 0 Deletion						
<input checked="" type="checkbox"/> <a href="#">1110002L01Rik</a>	Production...									
<a href="#">1110002O04Rik</a>	Production...									

Figure 2. Registering interest in a Micro-Injecting genes

2. Push the “Register Interest” button. That should create the plan. In this case the plan has been *Assigned*, because there were no other competing plans or production.

3. Edit the plan – if necessary – to specify what kind of allele you are intending to make. Clicking on the plan will bring up the plan edit page. This page has a section that allows the user to select an allele type. (Figure 3) This is very useful if you want to distinguish intent – for instance if you want to make it clear to other users that you intend to make a Knockout First Tm1a allele (or a Cre-Knockin). In the following I’ve checked the

'Knockout First Tm1a' box to make it clear for reporting purposes that I am going to make a Knockout First Tm1a (in case there are any deletion clones around).



▼ Allele Structure

*Bespoke Allele?*  
☐

*Knockout First Tm1a Allele?*  
☒

*Deletion Allele?*  
☐

*Cre Knock In Allele?*  
☐

*Cre Bac Allele?*  
☐

*Conditional Tm1c Allele?*  
☐

*Point Mutation Allele?*  
☐

*Conditional Point Mutation Allele?*  
☐

*Allele Symbol Override*


*Allele Comment*

Figure 3. Dialogue for changing details for an ES cell clone.

An ES Cell distribution centre (or a mouse clinic doing data entry on behalf of an ES-distribution centre) can indicate the progress of ES-QC for a plan by editing the data on the plan. This edit page is also accessible from the Plans module.

## 2b. Plans: Editing plans for Mouse Production.


The “Plans” module is selected by clicking the “Plans” tab . The resulting grid can be searched by Marker Symbol, Consortium, Production centre etc clicking on the column headers for these fields and choosing to “Filter” on that field. By default, the grid is filtered on the Users’s Production Centre (Figure 4).



Marker Symbol	Consortium	Production Centre	Status	Priority	Conditiona	Deletion	Cre Knock-	Cre BAC	Bespoke	Sub-Project
Atp6v0b	Sort Ascending Sort Descending Columns Filters	SI	Assigned	High	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Figure 4. Filtering plans in your consortium.

Clicking ‘Edit Form’ on a Plan row will produce the Plan edit-form for the selected plan. Entering a non-zero value in “# ES Cells starting QC” will cause the plan to change state from it’s previous state (“Assigned” or “Inspect ...” etc) to “ES QC in Progress” (Figure 5).



**ES Cell QC**

Number Of Es Cells Starting Qc:

Number Of Es Cells Passing Qc:

Es QC Comment:

**ES Cell QC**

Number Of Es Cells Starting Qc:

Number Of Es Cells Passing Qc:

Es QC Comment:

Marker Symbol: Cbx1

Consortium: MGP

Production Centre:

Status: Assigned - ES Cell QC In Progress

Phenotype Only: ☐

Mutagenesis Via Crispr/Cas9?: ☐

Marker Symbol: Cbx1

Consortium: MGP

Production Centre:

Status: Assigned - ES Cell QC Complete

Phenotype Only: ☐

Mutagenesis Via Crispr/Cas9?: ☐

Figure 5. Status change of ES clone to 1. “ES Cell QC in Progress” (above)  
2. “ES Cell QC Complete” (bottom)

Entering a non-zero value in “#ES Cells passing QC” will cause the plan to change state to “Assigned – ES QC Complete” (Figure 5). If the “# ES Cells passing QC” is set to zero, then the status of the Plan will change to “Assigned – ES QC Aborted”. This indicates that attempts to thaw / pass clones prior to Mouse production have failed for this consortium.


## 2c. Plans: Editing a plan to indicate crispr injection/ phenotyping

The Plans are not only used to show intent to produce mice from clones. They can also be used to show intent to:

- Create a mouse colony through mutagenesis via crispr injection (with or without a vector/Oligos).
- Phenotype an existing colony produced by another production centre.

These can be achieved through the edit plan page.

1. To show intent to create a mouse colony through mutagenesis via crispr injection tick the ‘Mutagenesis Via Crispr/Cas9’ tick box (figure 6) on the Plan and save.
2. To show intent to phenotype an existing colony produced by another production centre tick the ‘phenotype\_only’ tick box (figure 6) on the Plan and save



The screenshot shows a form with the following fields and values:

- Marker Symbol**: Cbx1
- Consortium**: MGP
- Production Centre**: WTSI (selected from a dropdown menu)
- Status**: Assigned - ES Cell QC In Progress
- Phenotype Only**: ☐ (unchecked)
- Mutagenesis Via Crispr/Cas9?**: ☐ (unchecked)

Figure 6. Showing intent to mutagenize Via Crispr injection or phenotype a mouse colony produced by another production centre

## 2d. Plans: Withdrawing interest in a gene, Inactivating a Plan, Deleting a Plan.

Consortia can

- Withdraw interest in a plan prior to the start of any mouse production
- Inactivate a plan with no *active* mouse production.
- Ignore Available Mice, hides mice from the public produced by this plan
- Delete an existing plan

By clicking on the obvious buttons on the plan's edit page (figure 7).

All of these changes will mean the plan does not 'compete' with other consortia looking to make plans on that gene.



The image shows a screenshot of a web form titled 'Other' with a dropdown arrow. It contains five checkboxes, each with a label above it: 'Report To Public' (checked), 'Ignore Available Mice' (unchecked), 'Recovery' (unchecked), 'Is Active' (checked), and 'Withdrawn' (unchecked). The form is enclosed in a light gray border.

Figure 7. Withdrawing interest in a gene, aborting, deleting and hiding produced mice from the public

### 3a. Attempts: Initiating Mouse Production for a gene.

iMits tracks Mouse production through Micro Injections of ES Cells or Mutagenesis Factors (via Crispr/Cas9).

NB. Before starting check the appropriate Plans and Targeting/HDR vector (for Mutagenesis Factor [see Adding Alleles/Vectors section]) have been created.

To record the start of a single Microinjection click the “Create MI Attempt” tab. This will give you two options.

- i. Select ES Cell or
- ii. Select/Create Mutagenesis factor

#### (i) Mouse Production using ES Cell.

ES Cells will be sent from the ES-distribution centre to the Mouse Clinic, where they will be prepared for Microinjection.

Click the ‘Select ES Cell’ button and search for clones that are visible for that gene, by either gene symbol or Clone name (Figure 8).

**Search for ES cells**

Search by marker symbol | Search by ES cell name

Enter marker symbol

Choose an ES cell clone to micro-inject

ES Cell	Marker Symbol	Pipeline	Mutation Subtype	LoxP Screen
EPD0738_6_B05	Atp6v0b	EUComm	conditional_ready	pass
EPD0738_6_C05	Atp6v0b	EUComm	targeted_non_conditional	not confirmed
EPD0738_6_C06	Atp6v0b	EUComm	targeted_non_conditional	not confirmed
EPD0738_6_D05	Atp6v0b	EUComm	targeted_non_conditional	not confirmed
EPD0738_6_D06	Atp6v0b	EUComm	targeted_non_conditional	no reads dete...
EPD0738_6_D08	Atp6v0b	EUComm	conditional_ready	pass

Figure 8. Searching for ES cell clones to choose for microinjection attempts

Choose the ES cell you have used for Microinjection (pay attention to the difference between conditional ready and targeted non-conditional clones, which is read from the IKMC targeting repository classification of these clones). I clicked on the row “EPD0738\_6\_B05”, which opens the Create form:



## (ii) Mouse Production using Mutagenesis Factor.

Click the 'Select/Create Mutagenesis Factor' button and type in the gene marker symbol (I entered Nxn) and click search. This will show you a list of exons which can be used to search for Crisprs/Crispr Pairs (Figure 9).

**Create mutation Factor**

Search for crisprs  
Enter marker symbol  
nxn Search

Exon	Sequence	Chr Start	Chr End
All	AAGACACGAAGACGATCTCCAGG	76398784	76398806
ENSMUSE00000392903	ACACGAAGACGATCTCCAGGCGG	76398787	76398809
ENSMUSE00000295885	ACCGTAGAAGGCGGCCAGGCTGG	76398887	76398909
ENSMUSE00000110352	ACGATCTCCAGGCGGTGCCGGG	76398795	76398817
ENSMUSE00000110354	AGACGATCTCCAGGCGGTGCCGG	76398793	76398815
ENSMUSE00000110358	AGCCCCGGCACGCGGGCGGGC...	76399074	76399096
ENSMUSE00000110357	ATCCTGGCCACCGCAGCCCGG	76399059	76399081
ENSMUSE00000110355	CCACCGCAGCCCGGCACGCGGG	76399067	76399089
ENSMUSE00000342546	CCACCGCAGCCCGGCACGCGGG	76399067	76399089
	CCACCGCCCGTCACGAGCTTGTC	76399008	76399030
	CCACCTCTCGCCACCGCCCGTC	76398997	76399019

**Mutagenesis Factor: Crispr Selection**

Sequence: Chr: Chr Start: Chr End: Add Crispr

Sequence	Chr	Chr Start	Chr End
ACGATCTCCAGGCGGTGCCGGG	11	76398795	76398817

Save Cancel

Figure 9. Crispr Selection using exon search

Alternatively you can manually enter (Figure 10) the 23 base sequence (this includes the PAM sites) and positional information of the Crispr.

**Create mutation Factor**

Search for crisprs  
Enter marker symbol  
nxn Search

Exon	Sequence	Chr Start	Chr End
All	AAGACACGAAGACGATCTCCAGG	76398784	76398806
ENSMUSE00000392903	ACACGAAGACGATCTCCAGGCGG	76398787	76398809
ENSMUSE00000295885	ACCGTAGAAGGCGGCCAGGCTGG	76398887	76398909
ENSMUSE00000110352	ACGATCTCCAGGCGGTGCCGGG	76398795	76398817
ENSMUSE00000110354	AGACGATCTCCAGGCGGTGCCGG	76398793	76398815
ENSMUSE00000110358	AGCCCCGGCACGCGGGCGGGC...	76399074	76399096
ENSMUSE00000110357	ATCCTGGCCACCGCAGCCCGG	76399059	76399081
ENSMUSE00000110355	CCACCGCAGCCCGGCACGCGGG	76399067	76399089
ENSMUSE00000342546	CCACCGCAGCCCGGCACGCGGG	76399067	76399089
	CCACCGCCCGTCACGAGCTTGTC	76399008	76399030
	CCACCTCTCGCCACCGCCCGTC	76398997	76399019

**Mutagenesis Factor: Crispr Selection**

Sequence: AGCCCCGGCACGCGGGCGGGC Chr: 11 Chr Start: 76399074 Chr End: 76399096 Add Crispr

Sequence	Chr	Chr Start	Chr End
ACGATCTCCAGGCGGTGCCGGG	11	76398795	76398817
AGCCCCGGCACGCGGGCGGGC...	11	76399074	76399096

Save Cancel

Figure 10. Manually entered crisprs

NB. At this point you can still edit your Crispr selection if incorrect.

Figure 11. Selected the nxn\_tv vector from the drop down

The MIDate is the microinjection date.

The Colony Name is a UNIQUE name CHOSEN BY YOU for the group of mice (the “mouse line”) that will arise from this microinjection. WTSI have colony names which are four-letter words such as “MECD”. UCD use colony names like “BL1253”. Typically this name is an “external” reference, determined by the name of the colony in the “real” mouse-tracking system available on the mouse clinic campus. This name will be transmitted to the outside world via portals, marts etc.

[illegible]

Figure 12. Adding microinjection details for a selected ES cell clone

Above I selected the UCD plan and entered the microinjection date and colony name. This successfully creates the Mi Attempt with status ‘Micro Injected in progress’ (figure. 13)

Micro-injection attempt created

ES Cell Details

Marker Symbol

Atp6v0b

ES Cell Name

EPD0738\_6\_B05

Figure 13. Green flag showing a microinjection attempt has been successfully entered

Create MI Attempt

Mouse Production

Genotype V

3b. Attempts: Indicating the production of chimeras/founders for a mouse line.

Mouse production attempts can be found and edited via the “Mouse Production” Grid. The grid can be searched by colony name, Gene symbol, or ES cell name. In Figure 14 I used colony name. Note also you can filter by production centre. The Grid can be directly filtered by consortium (by clicking on the “Consortium” column header and choosing to “Filter”). The Grid can be sorted by MI-date.

Here I’ve found my MI by filtering by colony:

Search for Micro-Injections

ES Cell Name, Gene Symbol, or Colony Name:

VIV1

Production Centre:

WTSI

Status:

Search

Clear

Micro-Injection Attempts

Everything

Transfer Details

Litter Details

Chimera Mating Details

QC Details

Summary

Edit In Form	Phenotype	Active Phenotype	Consortium	Production Centre	ES Cell	Marker Symbol	Allele symbol	MI Date	Status
<a href="#">Edit in Form</a>		0	BaSH	WTSI	EPD0738_6_B05	Atp6v0b	Atp6v0b <sup>tm1a(EU...</sup>	26-10-2012	Micro-injection in

Figure 14. Searching through microinjection attempts

You can edit this Mouse production attempt by clicking on the row, which will make the fields in the row “Editable”. Here I’ve clicked on the row and scrolled the view to the right, showing the field with “Total Blasts Injected” ready to type in:

Allele symbol	MI Date	Status	Colony Name	Distribution Centres	Blast Strain	Total Blasts Injected
Atp6v0b <sup>tm1a(EUCOMM)</sup>	26-10-2012	Micro-injection in progress	VIV1			0

Figure 15. Editing within the mouse production grid

You can also click on the “Edit in Form” link to the left, which will show the MI Attempt in a traditional HTML page.

‘Chimeras Obtained’ (ES Cell micro-injection) and ‘Founders Obtained’ (Mutagenesis Factor) statuses are equivalent for the different micro-injections. Clicking on the “Edit in Form” link will display slightly different forms depending on the different micro-injections. Figure16 shows the two forms and the required fields to change the Mi Attempt to “Chimeras Obtains” or “Founders Obtained”

### ES Cell

**Litter Details**  
Total Pups Born: 10  
Total Female Chimeras: 2  
Total Male Chimeras: 5  
Number of Males with 0-39% Chimerism: 3  
Number of Males with 40-79% Chimerism: 2

**Universal Details**  
Micro-Injection Date: 26/10/2012  
Status: Chimeras obtained  
Colony Name: VIV1  
Consortium: BaSH

### Mutagenesis Factor

**Crispr Injection Details**  
Number of Embryos Injected: 10  
Number of Embryos Survived: 10  
Number Transferred: 10

**Founder Litter Details**  
Number of Founder Pups Born: 2  
Number of Mutant Founders: 2  
Number of Founders Selected For Breeding: 2

**Universal Details**  
Micro-Injection Date: 26/10/2012  
Status: Chimeras obtained  
Colony Name: VIV1  
Consortium: BaSH

Figure 16. Changing status of mouse production to

Chimeras/Founder obtained

### 3c. Attempts: Indicating Genotype Confirmed mice for a mouse line.

A mouse line changes status to Genotype Confirmed when either of two fields is set to > 0: “Number of Chimeras with GLT from Genotyping” OR “Number of Het offspring”.

The following fields must also be filled in, so that details of the mouse line can be used to yield a correct strain name:

1. Test Cross Strain Name: the name of the initial test-cross strain used for breeding to establish chimerism.
2. Colony Background Stain Name: the name of the back-cross strain (eventually the colony background) for the genotype confirmed mouse colony.
3. Number of Chimeras with GLT from genotyping. The number of chimeras (ie parents) which had germ-line-transmitting pups, where the GLT was established by genotyping the pups. A non-zero number here will set the MI Attempt status to “Genotype Confirmed”.
4. Number of chimeras with x – y % chimerism. These bins are for information purposes only (they can be reported on).
5. Number of Het offspring. A non-zero number here will set the MI Attempt status to “Genotype Confirmed”.
6. Mouse Allele Type: If the ES Cell is a “tm1a” – ie knockoutfirst – then there is some chance that genotyping a mouse will show that the mouse has a “Tm1e” allele. This could be due to a number of reasons, but in this field users can record that the mouse was “tm1e”. If this field is left BLANK then the mouse allele type will be inherited from the allele of the ES Cell.

**Chimera Mating Details**

Colony Background Strain Name  
CS7BL/6N

Test Cross Strain Name  
CS7BL/6N

Date Chimeras Mated  
26/10/2012

Number Of Chimera Matings Attempted  
5

Number Of Chimera Matings Successful  
5

Number Of Chimeras With GLT From CCT  
5

Number Of Chimeras With GLT From Genotyping  
3

Number of Chimeras with 0-9% GLT  
0

Number of Chimeras with 10-49% GLT  
0

Number of Chimeras with 50-99% GLT  
5

Number of Chimeras with 100% GLT  
0

Total F1 Mice From Matings  
0

Number Of CCT Offspring  
0

Number Of HET Offspring  
7

Number Of Live GLT Offspring  
5

Mouse Allele Type  
a - Knockout-first - Reporter Tagged Insertion

Figure 17. Changing status of mouse production to Genotype Confirmed

### 3d. Attempts: Distribution Centres for Genotype Confirmed Mice.

A Genotype Confirmed line does not automatically inherit a distribution centre. **These must be added manually by the mouse production centre** (Figure 18).

**Distribution Centres**

Distribution centre	EMMA	Deposited material	Start date	End date	Action
<a href="#">Add Distribution Centre</a>					

Figure 18. Adding distribution centre to genotype confirmed mice

In Figure 18, I have clicked the “Add distribution centre” link, and chosen WTSI to act as an EMMA “node”, serving live mice. This tells the IMPC portal to refer users to the EMMA repository. It also tells the EMMA repository to expect to advertise this mouse line as order-able:

▼ *Distribution Centres*

Distribution centre	EMMA	Deposited material	Start date	End date	Action
<div>WTSI</div> <div></div>	<input checked="" type="checkbox"/>	<div>Live mice</div> <div></div>	<div></div>	<div></div>	<a href="#">remove</a>

[Add Distribution Centre](#)

Figure 15. Distribution Centre selection choices- WTSI as distributor under the EMMA framework.

IMPORTANT NOTE – to actually add the centre, you have to push the “UPDATE” button at the bottom of the web page:



The following distribution centre (see Figure 19) would be appropriate if WTSI sent the line to UCD for distribution, as part of the KOMP2 resource:

▼ *Distribution Centres*

Distribution centre	EMMA	Deposited material	Start date	End date	Action
<div>UCD</div> <div></div>	<input type="checkbox"/>	<div>Frozen embryos</div> <div></div>	<div></div>	<div></div>	<a href="#">remove</a>

[Add Distribution Centre](#)

Figure 19. Selecting KOMP repository for archiving of frozen embryos.

#### 4.A. Phenotyping Attempts: Registering a Genotype Confirmed Mouse for Modified Allele creation and Phenotyping

The Genotype Confirmed Atp6v0b mouse line (VIV1) can be registered for phenotyping by BaSH at WTSI, or by another consortium / centre.

Since it was made by BaSH / WTSI, to register it for BaSH at WTSI requires no further planning. Find the mouse line in the “Mouse Production” grid (see Figure 20) and click on the “Create” link in the “Phenotype” column:

Micro-Injection Attempts							
Everything	Transfer Details	Litter Details	Chimera Mating Details	QC Details	Summary		
Edit In Form	Phenotype	Active Phenotype	Consortium	<b>Production Centre</b>	ES Cell	Marker Symbol	Allele
<a href="#">Edit in Form</a>	<a href="#">Create</a>	0	BaSH	WTSI	EPD0738_6_B05	Atp6v0b	Atp6v0b

Figure 20. Selecting “Create” in the mouse production grid to start a phenotype attempt

You will be offered a Phenotype Attempt edit page (Figure 21), which allows you to Register a Phenotype Attempt for this mouse line.

ES Cell Details

Marker Symbol  
Atp6v0b

ES Cell Name  
EPD0738\_6\_B05

Universal Details

MI Attempt Colony Name  
VIV1

Colony Name  
VIV1-PA1

Consortium  
BaSH

Production Centre  
WTSI

Is Active  
☒

Distribution Centres

No Distribution Centres can be set until Status is Cre Excision Complete.

Rederivation


Rederivation Started  
☐

Figure 21. Creating a Phenotyping attempt, for the same Consortium and Production Centre as the generated mouse.



The COLONY NAME for the Phenotype Attempt is set BY YOU. **This name must match the colony name sent with your mouse clinic’s phenotyping data to the DCC.**

After I have pushed “create” for this Phenotype Attempt, it is in status “Registered”. You can find the Phenotype Attempt for update later by going to the “Phenotyping” tab:



You can search for an individual Phenotyping Attempts by entering their Gene in the search box (see Figure 22) . Alternatively, the columns of the grid will allow you to *filter* by consortium, production centre, colony name status etc. Below we have found the Phenotyping Attempt we just created for Atpv60b:

▼ Search for Phenotype Attempts

Gene Symbol:  
Atp6v0b

Production Centre: WTSI

Search

Clear

Phenotype attempts

Edit In Form	Colony Name	Consortium	Production Centre	Distribution Centres	Marker Symbol	Active?	Status	Rederivation started
<a href="#">Edit in Form</a>	VIV1-PA1	BaSH	WTSI		Atp6v0b	<input checked="" type="checkbox"/>	Phenotype Attempt Registered	<input type="checkbox"/>

Figure 22. Searching for Phenotyping Attempts already created

Clicking on this row brings up the Phenotype Attempt edit form, which allows updates to the status of the Phenotype Attempt.

You can view the Production Graph (Figure 23) of this whole gene by selecting it in the “Gene Selection” tab, and clicking on the “Production Graph” link. You will see the gene, the plan, the mouse and the phenotype attempt all linked up:




Figure 23. The Production Graph for the mouse and phenotype attempt, visible via the Gene Selection Page.

Phenotype Attempts: Having multiple Phenotype Attempts for a single mouse line (from different consortia or phenotyping centres)

What if BCM also wanted to Phenotype this same mouse line (VIV1)? Well, they could! But it would require the creation of a new PLAN.

Recall that the Plan we created for making the mouse “VIV1” was for BaSH / WTSI. First we would have to create a NEW plan – this time for phenotyping – for BaSH / BCM: I have created a new PLAN for Atp6v0b for BaSH / BCM using the gene selection tab (see Figure 24).

Notice that it is in status “Inspect – GLT Mouse” because iMITS has observed the existing GLT mouse at WTSI, and is questioning BCM’s intention to produce. Note also that iMITS is showing the existing GLT mouse and the Phenotype Attempt already registered at WTSI.

The screenshot shows the 'Filter Genes' interface. At the top, there is a search bar with 'Atp6v0b' entered. Below it, there are two dropdown menus: 'MIs for Consortium:' and 'MIs at Production Centre:'. The 'MIs for Consortium:' dropdown is set to 'BaSH' and the 'MIs at Production Centre:' dropdown is set to 'BCM'. Below these, there is a table titled 'Please Select the Genes You Would Like to Register Interest In\*'. The table has columns for 'Gene', 'Production History', '# IKMC Projects', '# Clones', 'Non-Assigned Plans', 'Assigned Plans', 'Aborted MIs', 'MIs in Progress', 'Genotype Confirmed', and 'Phenotype Attempt'. The first row shows 'Atp6v0b' with a 'Production Graph' link, 4 IKMC Projects, 2 Conditional and 9 Targeted Traps, and a link to '[BaSH:BCM:Inspect - GLT Mouse]'. The 'Phenotype Attempt' column shows '[BaSH:WTSI:1]' and '[BaSH:WTSI:1]'.

Figure 24. Creating a new plan for BaSH / BCM, for the Atp6v0b gene

Now that we have made a plan for BaSH / BCM, the NEXT thing to do is to actually FIND the mouse and make a Phenotype Attempt attached to the BaSH / BCM plan. Go to the Mouse Production Tab, and click on the ‘Create’ link in the Phenotype column (Figure 25).

Micro-Injection Attempts							
Everything	Transfer Details	Litter Details	Chimera Mating Details	QC Details	Summary		
Edit In Form	Phenotype	Active Phenotype	Consortium	Production Centre	ES Cell	Marker Symbol	Allele
<a href="#">Edit in Form</a>	<a href="#">Create</a>	0	BaSH	WTSI	EPD0738_6_B05	Atp6v0b	Atp6v0b

Figure 25. Finding the original Atp6v0b mouse prior in order to create a BaSH / BCM phenotyping attempt.

THIS time you want to create a Phenotype Attempt for consortium: BaSH, production centre: BCM (Figure 26) –

ES Cell Details

Marker Symbol

Atp6v0b

ES Cell Name

EPD0738\_6\_B05

Universal Details

MI Attempt Colony Name

VIV1

Colony Name

VIV1-PAforBaSH-BCM

Consortium

BaSH

Production Centre

BCM

Is Active

☒

Figure 26. Selecting BaSH / BCM when creating a Phenotyping attempt for a BaSH/WTSl mouse.

After you create the PA (push “Create” at the bottom of the screen) then you can check your work by looking at the Production Graph for this gene via the row for the gene in GeneDetails page (Figure 27):





Figure 27. The Production graph for Atp6v0b after creation of the BaSH/BCM Phenotyping attempt on the BaSH/WTSl mouse.

You will see that there are two plans in the system: one for BaSH/WTSl which resulted in a Genotype Confirmed tm1a mouse (and a phenotype attempt registered), and the other for BaSH/BCM, which resulted in a *separate* Phenotype Attempt for the *same* mouse.


#### 4b. Phenotyping Attempts: Indicating the start / finish of rederivation

Sometimes a mouse clinic will have to rederive a colony after the reception of cryo-preserved stock (embryos, sperm). iMITS allows the optional indication of that process by the checking of the “Rederivation Started” and “Rederivation Complete” checkboxes on the Phenotype Attempt edit page.

To get to this page, first navigate to the Phenotyping Attempts tab:



and search for the attempts we've just created by entering the gene name (remove the 'Production Centre' filter (Figure 28):

A screenshot of a web application interface. At the top is a search form titled 'Search for Phenotype Attempts'. It has a text input field for 'Gene Symbol:' containing 'Atp6v0b' and a dropdown menu for 'Production Centre:'. Below the input field are 'Search' and 'Clear' buttons. Below the search form is a table titled 'Phenotype attempts'. The table has columns: 'Edit In Form', 'Colony Name', 'Consortium', 'Production Centre', 'Distribution Centres', 'Marker Symbol', and 'Active?'. There are two rows of data. The first row has 'Edit in Form' as a link, 'VIV1-PA1' as the colony name, 'BaSH' as the consortium, 'WTSI' as the production centre, an empty distribution centres field, 'Atp6v0b' as the marker symbol, and a checked 'Active?' checkbox. The second row has 'Edit in Form' as a link, 'VIV1-PAforBaS...' as the colony name, 'BaSH' as the consortium, 'BCM' as the production centre, an empty distribution centres field, 'Atp6v0b' as the marker symbol, and a checked 'Active?' checkbox.

Edit In Form	Colony Name	Consortium	Production Centre	Distribution Centres	Marker Symbol	Active?
<a href="#">Edit in Form</a>	VIV1-PA1	BaSH	WTSI		Atp6v0b	<input checked="" type="checkbox"/>
<a href="#">Edit in Form</a>	VIV1-PAforBaS...	BaSH	BCM		Atp6v0b	<input checked="" type="checkbox"/>

Figure 28. Finding the Phenotyping Attempts for Atp6v0b in the Phenotyping grid, in order to update the BaSH/BCM attempt.

Notice that we're now finding both the Phenotype Attempts we created: one at WTSI and one at BCM. Click on the “*Edit in Form*” link for the BCM row. This will bring up the Phenotype Attempt Edit form. The middle of the form (Figure 29) has checkboxes for “Rederivation Started” and “Rederivation Complete”:



**Rederivation**

Rederivation Started ☒

Rederivation Complete ☒

**Cre Matings**

Cre-deleter Strain

MGI:3046308: Hprt<tm1(CMV-cre)Brd>

Number Of Cre Matings Successful

0

Mouse Allele Type

Figure 31. Specifying the cre-deleter strain to alter the status of the Phenotyping attempt to “Cre Excision Started”.

### Phenotyping Attempts: Indicating the finish of Cre Excision

The *finish* of Cre-Excision happens when a mouse clinic has genotyped the offspring of the cre-mating and determined which matings have successfully removed the promoter / floxed the critical exon(s).

To change to this status, the mouse clinic must indicate *all* of:

- A positive number of “cre matings successful” (this is a relic and will be replaced with a checkbox)
- A Mouse Allele Type (this should be a “b” for a conditional allele and a “.1” for a deletion)
- A Colony Background Strain. This is the colony background of the stain which will be phenotyped.

We have made choices in the form shown below (Figure 32): pushing “Update” will change the status of the Phenotype Attempt to “Cre Excision Complete”.



### Phenotyping Attempts: Distribution Centres for Cre-Excised Mice.

A Cre-Excised mouse line. can be distributed as a resource available to the scientific community, just like the original conditional – ready mouse line.

Check the “Distribution Centres” panel for your Phenotype Attempt: you can add a distribution centre there - and indicate if it is to be distributed via the EMMA network - in the same way as the original mouse lines (Figure 34).

Distribution centre	EMMA	Deposited material	Start date	End date	Action
BCM	<input type="checkbox"/>	Frozen embryos			<a href="#">remove</a>

[Add Distribution Centre](#)


Figure 34. Specifying / adjusting the Distribution Centre for the Cre-Excised mouse on the Phenotype Attempt edit page.

### Phenotyping Attempts: Indicating the start and end of Phenotype Data Transfer

The iMITS status currently labeled “Phenotyping Started” and “Phenotyping Complete” represent the start and finish of phenotyping data transfer between the mouse clinic and the Pheno-DCC - – the body responsible for the compilation and quality-control of phenotype data from the entire KOMP2 / IMPC effort.

Mouse clinics cannot directly indicate the start and end of data transfer inside iMITS: This signal must be sent to iMITS directly by the data wranglers at the Pheno-DCC (Figure 35). Currently this coordination requires that the Phenotype Attempt Colony name recorded in iMITS match the Colony Name attached to the XML containing the procedure / parameter results in the data transferred to the Pheno-DCC:





1. Data packet from mouse clinic contains PA colony name
2. Pheno DCC reads data packet
3. Pheno DCC updates iMITS to change status to "Phenotype Data Transfer Started" (and "Complete")

Figure 35. Coordination between iMITS, Data coming from the Mouse Clinic and the PhenoDCC to indicate the flow of data via "Phenotyping Started / Finished" statuses

## 4. Adding Alleles/Vectors

The IKMC Targeting Repository (Targ Rep) Stores all the mutant ES Cells and Targeting Vectors made by the pipelines in the International Knockout Mouse Consortium (IKMC). The Targ Rep has also been adapted to stores the Targeting Vectors and Oligos created for Crispr/Cas9 microinjections.

Here we will cover how to add these new Targeting Vectors and Oligos.

**Targ Rep will infer most of the alleles structure from the information provided when entering the Targeting Vectors and Oligos. Therefore it is IMPORTANT that the information entered is correct to prevent an incorrect allele from being inferred.**

Targ Rep's navigation bar has a Crispr and Cell section. The Crispr section has two options 'Targeted Alleles' and 'HDR Alleles'. Click on 'Targeted Alleles' to add Vectors and click on 'HDR Alleles' to add Oligos.

Crisprs Targeted Alleles HDR Alleles

### Adding Targeting Vectors (Targeted Alleles)

Click on 'Targeted Alleles' and then select 'new'. The following form will pop up.

Here I have filled in the gene as Nxn, the genomic information, design details, the targeting vectors Genbank file and the name of the targeting vector (figure 36). **NB.** You can enter more than one vector name if you have multiple copies of the vector. It is important that the Genbank file is accurately annotated as this will be used to infer the allele's structure.

*Create an allele*

Gene Filter

Nxn

Genomic Position

Chromosome

11

Strand

+

Assembly

GRCm38

Design Details

Mutation Method

Targeted Mutation

Mutation Type

Conditional Ready

Mutation Subtype

Domain Disruption

Backbone

L3L4\_pD223\_DTA\_T\_spec

Cassette

pL1L2\_GT2\_LF2A\_nEGFP0\_T2A\_CreERT\_puro\_dr

Cassette type

Promotorless

Genbank files Details

Targeting Vectors Genbank File

```
LOCUS      KO-  
first_condition_ready_114372_MGI:2677836  
21996 bp    dna      circular UNK  
DEFINITION Mus musculus targeted KO-first,  
conditional ready, lacZ-tagged  
mutant vector Dolk targeting  
project(s): 114372  
ACCESSION   unknown  
SOURCE      Mus musculus  
ORGANISM    Mus musculus  
Eukaryota; Metazoa; Chordata;  
Craniata; Vertebrata; Euteleostomi;  
Mammalia; Eutheria;  
Euarchontoglires; Glires; Rodentia;  
Sciurognathi; Muroidea; Muridae;  
Murinae; Mus.  
COMMENT     cassette:
```

*Targeting Vectors*

Targeting Vector	Intermediate Vector
nxn_tv	

[Add a targeting vector](#)

Figure 36.  
Creation of a  
targeting  
vector

Click 'save'. Now this vector will appear when recording the microinjection of a Mutation Factor (Crisprs and vector)

## HDR Allele (Oligos)

Click on 'HDR Alleles' and then select 'new'. The following form will pop up.

Here I have filled in the gene as Nxn and the genomic information. The design details default to 'Point Mutation' and I have entered the sequence of the oligos and the name of the vector (figure 37). **NB.** You can enter more than one vector name if you have multiple copies of the same vector. It is important that the Oligos sequence is accurately entered as this will be used to infer the allele's structure.

Gene Filter

Nxn

Genomic Position

Chromosome

11 ▼

Strand

+ ▼

Assembly

GRCm38

Design Details

Mutation Method

Targeted Mutation ▼

Mutation Type

Point Mutation ▼

Mutation Subtype

Point Mutation ▼

Subtype Description

Backbone

Oligos Details

Sequence

ccactggccgctcggttttacattaattaacaacttttctatacaagtt  
gacagatctatagtgtcacctaaatccaaaaaacggccaacattg  
gcccggttttttccggtttatctgtttaactcggccgctctagcctc  
gaggctagaactagtggatctcgagccccagctgggttctttccgcctc  
agaagccatagagcccaccgcatccccagcatgctgctattgtcttc

Targeting Vectors

Targeting Vector	Report to Public	
nxn_hdr_oligos	<input type="checkbox"/>	<a href="#">Remove</a>

[Add a targeting vector](#)

Figure 37.  
Creation of  
Oligos

Click 'save'. Now this oligos will appear when recording the microinjection of a Mutation Factor (Crisprs and oligos)

## 5. GLOSSARY

### IKMC

The members of the International Knockout Mouse Consortium (IKMC) have worked together to mutate all protein-coding genes in the mouse using a combination of gene trapping and gene targeting in C57BL/6 mouse embryonic stem (ES) cells. The IKMC includes the following programs: The Knockout Mouse Project, The European Conditional Mouse Mutagenesis Program, The North American Conditional Mouse Mutagenesis Project, The Texas A&M Institute for Genomic Medicine, and the EUCOMMTools Program. See [www.knockoutmouse.org](http://www.knockoutmouse.org) for details of IKMC products, alleles and contributors.

### Genes

A mouse has (about) 20,000 *genes* inside every cell in its body, each encoded by a stretch of DNA. Each gene is labelled by an *MGI Accession Id*, which can be thought of as the logical key of the gene. The *Marker Symbol* is a human-readable string that labels the gene.

### ES Cells

iMITS contains a list of all IKMC ES Cells for a particular gene. When a user production center wants to start mouse creation, they must first pick a gene, and then an ES Cell from a list of possible ES Cells for that gene.

### Consortia

Mouse Production is funded by various Government and EU-related bodies known as Consortia. Each Consortium is granted money and held accountable for certain mouse-production goals. Examples of consortia include BaSH, the MGP, Phenomin, Helmholtz-GMC etc.

### Production Centres

Mouse production occurs in physical locations called Production Centres. These correspond to actual labs that receive ES Cells (from the ES Cell distribution centres) and then proceed to check the ES Cells and breed mice derived from them. Examples of production centres include WTSI, BCM, MRC Harwell etc.

One Consortium can have work done by many mouse production centres, and one mouse production centre can do work on behalf of many different consortia. The connection between the two is held by the MI Plans – see below.

## MI Plans (or just Plans, really)

The MI Plan encodes the *intention* of a consortium to produce a mouse for a particular gene at a particular production centre. The nature of the mutant mouse can be carefully specified by flags on the MI Plan.

## MI Attempts

These records are linked to MI Plans and ES Cells, and represent the creation / progress of colonies (groups of genetically identical) mice based on a single mutant ES cells. The records are linked to a single ES Cell, and proceed through statuses:

## Mutagenesis Factors

These are a combination of Crisprs/Cas9 and vectors that will be micro-injected into a mouse embryo to induce a mutation in the targeted gene.

## Phenotype Attempts (better called “Modified Alleles”)

These records are linked to MI Attempts (and independently to MI Plans if necessary) and represent the import of an existing mutant mouse, the modification (via cre-excision) of the allele in that mouse, and the start and end of phenotype data transfer to the DCC for that mouse.

## Statutes in iMITS

The overall progress of mouse production from planning to end of phenotyping is captured by status on Plans, MIAttempts and Phenotype Attempts for a gene:

Entity	Status	Description
Plan	Assigned	There are no other assigned plans for this gene (no other Consortia and Centres have plans for this gene)
	Assigned – ES QC In Progress	ES Cells for this gene have started QC, either at the ES-cell-distribution centre or at the mouse clinic.  Triggered by selecting non-zero ES Cells entering QC on this gene.  Note – simply entering ES Cells for QC will cause the plan to enter this state, even if it was not assigned before.
	Assigned – ES QC Complete	ES Cells for this gene have (successfully) finished QC prior to mouse production.  Triggered by selecting non-zero ES Cells finished QC on this gene.
	Assigned – ES QC Aborted	All ES Cells starting QC for this gene have failed to pass QC. The plan is <i>not yet</i> inactive or withdrawn.
	Inspect – Conflict	There are other plans from other consortia for this gene, but no visible mouse production. Note – the MI Plan can still be changed to status “Assigned” by starting ES QC, or by starting mouse production (MI)
	Inspect – MI Attempt	There are other visible (active) mouse production attempts for this gene. Note – the MI Plan can still be changed to status “Assigned” by starting ES QC, or by starting mouse production (MI)
	Inspect – MI Attempt / GLT Mouse	There are other visible Genotype Confirmed Mice for this gene. Note – the MI Plan can still be changed to status “Assigned” by starting ES QC, or by starting mouse production (MI)
MIAttempt	MI In Progress	Mouse production has commenced with Microinjection for this gene. There is a tracked mouse ‘colony’ – group of mice.
	Chimeras Obtained	The mouse colony has produced Chimeras. Triggered by entering > 0 male Chimeras for the MI Attempt.
	Founder Obtained	The mouse colony has produced founders from

		the micro-injection of the Mutagenesis Factor. Triggered by entering >0 Number of Mutant Founders
	Genotype Confirmed (ES Cell)	The mouse colony has produced F1 pups which have been assayed and found genotype confirmed. All centres except WTSI: Triggered by > 0 numbers in either of these fields: "Chimeras with glt offspring established by genotype confirmation" or "Number of het offspring".
	Genotype Confirmed (Mutagenesis Factor)	Genotype confirmation of the F1 (bread from the founders). CURRENTLY NOT POSSIBLE
	MI Aborted	
Phenotype Attempt	Phenotype Attempt Registered	This is the initial status on creation of a PhenotypeAttempt as a child of an MI Attempt
	Rederivation Stated / Finished	Triggered by the selection of a checkbox on the Phenotype Attempt.
	Cre Excision Started / Complete	Cre-started – triggered by selection of >0 for "Number of Cre Matings". Cre-complete: triggered by choice of "Mouse Allele Type" to "b".
	Phenotype Data Capture Started / Complete	These are only set-able by the DCC – that is, they cannot be set by the Mouse Clinic